

customer:

date: 29.10.2024

Harman Professional Denmark ApS.  
Olof Palmes Allé 44

8200 Aarhus N  
Denmark

project No.: 24000108-5

inspection order: colorimetric measurement according DIN EN 13032-4:2015

- evaluation of colour temperature
- evaluation of colour rendering

measurement procedure: - measurement with doublemonochromator 300-800nm

test sample specification:

luminaire: LED fixture for stage and studio, direct beam

description: MAC Viper XIP - medium beam

Serial No.: 15077051976

DMX settings:

dimmer:	65535
zoom:	22528
focus:	30080
color wheel:	1 (CRI-Filter)
fans:	full fan


test conditions:

stabilisation time: 30 min.  
 operating voltage: 230 V AC (const. <0.1%)  
 ambient temperature: 25 ° C ± 1°C


colorimetric data:

Color Temp. [K]	x- Coordinate	y- Coordinate	$\Delta uv$	CRI / Ra
5783	0,3273	0,3089	-0,0151	87,2

TM 30 Rf	TM 30 Rg	TLCI	CQS
80	108	68	82

colour temperature at CIE - colour triangle:

spectral light distribution:



The company ILUmetriX insures that all statements are made by the best of knowledge, according to European standards.

Meschede, 29.10.2024

*F. Siepe*

F. Siepe (Management Board)



annex: page 4 test sample  
page 5 test equipment

testsample:

dimensions [L x W x H]: 430 x 270 x 760 mm



*Test Equipment Spectrometry - ILUmetriX GmbH*

**IDR 300-PSL Spectroradiometer System, Bentham Instruments Ltd. UK**

Monochromator configuration: Symmetric, Single Czerny-Turner  
Monochromator focal length: Single 300mm, double 600mm  
Bandwidth: Software controlled motorised slit  
Number of gratings: 1-3 mounted on turret

Resolution:		<b>Single</b>	<b>Double</b>
	2400g/mm	-	0,075 nm
	1200g/mm	0,3 nm	0,15 nm
	400g/mm	0,9 nm	0,45 nm
Dispersion:	2400g/mm	-	0,68 nm/mm
	1200g/mm	2,7 nm/mm	1,35 nm/mm
	400g/mm	8,1 nm/mm	4,05 nm/mm
Wavelength accuracy:	2400g/mm	-	±0,1 nm
	1200g/mm	±0,2 nm	±0,2 nm
	400g/mm	±0,6 nm	±0,6 Nm

**Detection Electronics**

Current amplifier: Six decade trans-impedance amplifier  
Gain ranges:  $10^{10} - 10^5$  V/A  
Analogue digital converter: 100ms Integration

**Input Optics**

Light Transport: Flexible quartz fibre bundle, 200 - 1400 nm  
Diffuser cosine response:  $f^2$  error <1% 200 - 1100 nm  
Sphere cosine response:  $f^2$  error <5% 200 - 3000 nm  
Telescope, fields of view: 1.7 + 11 mrad  
Tel- measurement range: 200mm to 50m

**Calibration Standards**

Irradiance: CL6 Irradiance Standard  
Serial No: 13449/1  
Radiance: SRS12 Radiance Standard  
Serial No: 13559/1  
Traceability: NPL, UK

**BTS256-EF Lux- and Luminous Color Lightmeter, Gigahertz-Optik GmbH**

Serial No.: 50491